



2021 PRODUCT GUIDE







MPK

MDK Co., Ltd. is a manufacturer of precision dental handpieces and accessories in South Korea. MDK has consistently elevated its manufacturing capability since the foundation and is now ready to be recognized as one of the leading producers of top-quality but competitively priced dental handpieces and accessories.

MDK is now looking towards expanding its distribution network and increasing its brand value through a global strategy of progressively developing new handpiece series, manufacturing more value-added products, providing personalized services, and opening new international markets and industries.

MDK products are KGMP,CE and FDA 510(K) approved.

For information on MDK handpieces please visit www.mdk-dental.com.

Precision Crafted Performance

For your safe and effective preparations



Product Features

Crafted Titanium Body

MDK is specialized in producing dental handpieces in solid titanium body. Titanium is about 30 % lighter than conventional materials for dental handpieces, also strong against repeated sterilization.

The light weight reduces clinicians wrist fatigue significantly and non-slip surface of solid titanium makes the whole procedure safer.

Ceramic Ball Bearing

MDK only use Ceramic ball bearings which generates less wear and less temperatures. Ceramic balls increase the bearing's speed capabilities and life span compared to the Steel ball bearings.

Advanced Spindle-producing Technology

Bur eccentricity causes bur fly-out, bur wobbles and failure of turbine. MDK's spindle(where the bur inserted) is made of Carbide Steel which is very hard metal and the chuck extracting force is ≥30N.

Pressure Control System

MDK's pressure control system ensures only stable amount of air-pressure is provided while operation so the bearings would not be damaged.

Multi Water Spray

MDK's turbines offer from single to multi-port water spray that provide ideal cooling effect for the entire operation.

Fiber Optic System

MDK's Fiber Optic system enhances operator visibility.

Fiber optics are made from solid glass rod instead of a collection of optic fibers.

MDK's LED light bulb generates a brighter and whiter light.

Noise-Free Operation

MDK's advanced balancing technology by high-end balancing equipment provides very quiet operation

And noise level is less than 60dB(A)

Scratch Resistant treatment

MDK 's surface treatment technology offers scratch resistant and comfortable grip.

I TITANIUM

PEDIATRIC HEAD

STANDARD HEAD

TORQUE HEAD

45° HEAD



Ø11.2

14



414 4Ø11.2





Model: MOKO 560N LED Order code: 560NLS Cartridge code: N600



Model: MOKO 760N LED Order code: 760NLT Cartridge code: K8000



Model: MOKO TN45N LED Order code: TN45NL Cartridge code: N600

NSK®

COMPATIBLE MDK HANDPIECES





Model: MDK-2C-N1 Order code: MDK-N1 FIBER OPTIC 6-PIN TYPE





Model: MDK-2C-N4 Order code: MDK-N4 E- LIGHT GENERATOR TYPE

I NON-TITANIUM

STANDARD HEAD



Model: MOKO 520N LED Order code: 520NLS Cartridge code: N600

TORQUE HEAD



Model: MOKO 720N LED Order code: 720NLT Cartridge code: K8000

45° HEAD



Model: MOKO N45N LED Order code: N45NL Cartridge code: N600

KAVO®

COMPATIBLE MDK HANDPIECES





Model: MDK-2C-N1 Order code: MDK-N1 FIBER OPTIC 6-PIN TYPE





Model: MDK-2C-N4 Order code: MDK-N4 E- LIGHT GENERATOR TYPE

NSK is registered trademarks of NAKANISHI INC.

PEDIATRIC HEAD

STANDARD HEAD

TORQUE HEAD

Ø12.5

45° HEAD



Model: MOKO 760K LED PEDIATRIC Order code: 760KLP Cartridge code: KPEDI



Model: MOKO 560K LED Order code: 560KLS Cartridge code: N600



Model: MOKO 760K LED Order code: 760KLT Cartridge code: K8000



Model: MOKO TN45K LED Order code: TN45KL Cartridge code: N600

KAVO®

COMPATIBLE MDK HANDPIECES





Model: MDK-2C-K1 Order code: MDK-K1 FIBER OPTIC 6-PIN TYPE





Model: MDK-2C-K4 Order code: MDK-K4 E-LIGHT GENERATOR TYPE

I NON-TITANIUM

STANDARD HEAD



Model: MOKO 520K LED Order code: 520KLS Cartridge code: N600

TORQUE HEAD



Model: MOKO 720K LED Order code: 720KLT Cartridge code: K8000

45° HEAD



Model: MOKO N45K LED Order code: N45KL Cartridge code: N600

KAVO®

COMPATIBLE MDK COUPLERS





Model: MDK-2C-K1 Order code: MDK-K1 FIBER OPTIC 6-PIN TYPE





Model: MDK-2C-K4 Order code: MDK-K4 E-LIGHT GENERATOR TYPE




| NON-OPTIC with Coupler

KAVO® TYPE

TORQUE HEAD

- Model: MOKO 710K
- Order code: 710K Cartridge code: K8000



Model: MDK-2C-K01/02Order code: MDK-K01/02





Model	Cartridge code	Head type
*TITANIUM MOKO 750K	K8000	Torque head
*TITANIUM MOKO 550K	N600	Standard head
MOKO 710K	K8000	Torque head
MOKO 510K	N600	Standard head

INON-OPTIC with Coupler

NSK® TYPE

TORQUE HEAD

- Model: MOKO 710N
- Order code: 710N Cartridge code: K8000







- Model: MOKO K45N
- Order code: K45N Cartridge code: K8000



Model: MDK-2C-N01Order code: MDK-N01



Model	Cartridge code	Head type
*TITANIUM MOKO 750N	K8000	Torque head
*TITANIUM MOKO 550N	N600	Standard head
MOKO 710N	K8000	Torque head
MOKO 510N	N600	Standard head
MOKO K45N	K8000	45° head

SIRONA® BIEN-AIR®

COMPATIBLE MDK HANDPIECES

SIRONA® TYPE

STANDARD HEAD TORQUE HEAD



Model: MOKO 560S LED Order code: 560SLS Cartridge code: N600

Model: MOKO 760S LED Order code: 760SLT Cartridge code: K8000

BIEN-AIR® TYPE

STANDARD HEAD TORQUE HEAD



Model: MOKO 560B LED Order code: 560BLS Cartridge code: N600 Model: MOKO 760B LED Order code: 760BLT Cartridge code: K8000

SIRONA® BIEN-AIR®

COMPATIBLE MDK COUPLERS





Model: MDK-2C-S2 Order code: MDK-S2 FIBER OPTIC 6-PIN TYPE





Model: MDK-2C-B2 Order code: MDK-B2 FIBER OPTIC 6-PIN TYPE

Specifications are subject to change without notice

BORDEN® MIDWEST®

COMPATIBLE MDK HANDPIECES WITH PRESSURE CONTROL SYSTEM *





STANDARD/MULTI

TORQUE HEAD

STANDARD/SINGLE



Model: MOKO 230BT Model: MOKO 230B Order code: 230BT Order code: 230B Cartridge code: K8000

Model: MOKO 220BT Order code: 220BTS Cartridge code: N600 Model: MOKO 220B Order code: 220BS Cartridge code: N600

Cartridge code: K8000

Model: MOKO 250B Order code: 250B Cartridge code: MPANA

* WHAT IS PRESSURE CONTROL SYSTEM?

Sometimes dental clinics use higher air pressure than recommended level which leads to early bearing failure. Pressure control system ensures only stable amount of air-pressure is provided while operation so the bearings would not be damaged.

BORDEN® MIDWEST®

COMPATIBLE MDK HANDPIECES WITH PRESSURE CONTROL SYSTEM *







Model: MOKO 420MT Order code: 420MT Cartridge code: N600

Model: MOKO 420M Order code: 420M Cartridge code: N600

TORQUE HEAD



Model: MOKO 430MT Order code: 430MT Cartridge code: K8000

Model: MOKO 430M Order code: 430M Cartridge code: K8000

STANDARD/SINGLE



Model: MOKO 450M Order code: 450M Cartridge code: MPANA

* WHAT IS PRESSURE CONTROL SYSTEM?

Sometimes dental clinics use higher air pressure than recommended level which leads to early bearing failure. Pressure control system ensures only stable amount of air-pressure is provided while operation so the bearings would not be damaged.

BORDEN® MIDWEST®

COMPATIBLE MDK HANDPIECES WITH SELF-LIGHT GENERATOR



MDK's GT handpieces are produced by our precision engineering and show high performance and long-term reliability. GT handpieces are produced with a built-in power generator that delivers quality LED light when connected to any standard non-optic air tubing.

FEATURES OF MDK DENTAL CARTRIDGES

MDK OFFERS TOP QUALITY DENTAL REPLACEMENT PARTS FOR A WIDE RANGE OF THE MAJOR BRANDS.OEM PARTS CAN BE PRODUCED UPON REQUEST.

DIMENSIONAL ACCURACY: \leq 0.006MM BUR RUNOUT ACCURACY: \leq 0.006MM ACCURACY OF THE BUR GUIDE: 0.002MM CHUCK EXTRACTING FORCE: \leq 30N PUSH BUTTON OPENING FORCE: \leq 45N

HIGH PERFORMANCE CERAMIC BEARINGS
DURABLE CARBIDE MADE BUR GUIDE
ACCURATELY BALANCED BY HIGH-END BALANCING EQUIPMENT



HEAD SIZE	ORDER CODE	COMPATIBILITY
MDK TORQUE HEAD	K8000	KAVO® 8000B
MDK STANDARD HEAD	N600	NSK® TI-MAX 600
MDK PEDIATRIC HEAD	KPEDI	N/A
MDK STANDARD HEAD	MPANA	NSK® PANA MAX II

MDK LOW SPEED HANDPIECES

1:1 Constant speed contra-angle



Model: CA-PL1 Order code: FC001 Max. Motor Speed: 40,000 RPM Fiber optic/For CA burs

20:1 Reduction speed contra-angles for implantology and surgical applications



Model: CA-22PM Order code: FC004 Non-optic



Model: CA-PL20 Order code: FC004-S Non-optic



Model: CA-PL20X Order code: FC003 Fiber optic



Model: CA-PL30 Order code: FC003-S Fiber optic

1:5 Increasing speed contra-angle



Model: FG-PL15 Order code: FC005 Max. Motor Speed: 40,000 RPM Fiber optic/For FG burs

FG-PL15's Long service life, quiet operation with the highest quality bearings and gears allow excellent treatment results.

- Push-button chuck
- · Can be connected to all ISO motors
- Glass rod for constant light quality and visibility

MDK LOW SPEED HANDPIECES

1:1 Non optic contra angles

EXTERNAL SPRAY



Model: CA-20L Order code: C-001 LATCH TYPE MAX. SPEED: 25,000 RPM



Model: CA-20P Order code: C-002 PUSH BUTTON TYPE MAX. SPEED: 25,000 RPM



Model: CA-20PB Order code: C-003S PUSH BUTTON TYPE MAX. SPEED: 25,000 RPM

INTERNAL SPRAY



Model: CA-20Pi Order code: C-003 PUSH BUTTON TYPE MAX. SPEED: 25,000 RPM

1:1 Non optic straight handpieces

EXTERNAL SPRAY



Model: MD-4 Order code: M-001 MAX. SPEED : 25,000 RPM



Model: MD-5 Order code: M-003S MAX. SPEED: 25,000 RPM

INTERNAL SPRAY



Model: MD-4i Order code: M-002 SMOOTH GRIP MAX. SPEED: 25,000 RPM



Model: MD-5i Order code: M-004 STRONG GRIP MAX. SPEED : 25,000 RPM

MDK LOW SPEED E-TYPE AIR MOTORS

Fiber optic air motor with Internal spray (6-pin type connection)





Model: T260i-LED Order code: AM1001

Non optic air motor with Internal spray (2 and 4-hole type)





Model: 260i-B2 Order code: AM1002





Model: 260i-M4 Order code: AM1003

Non optic air motor with external spray (2 and 4-hole type)







Model: 260E-B2 Order code: AM1004






Model:260E-M4 Ordercode:AM1005







604, Ilsan-techno-town, Ilsan-ro 138, Gyeonggi-do, 10442 South Korea

E-MAIL: mdk@mdk-dental.com WEBSITE: WWW.MDK-DENTAL.COM

TEL: +82-(0)31-906-5601~2 FAX: +82-(0)31-5600

